

T-CLASSIC 9

9-Liter Portable Autoclave

Fast. Compact. Portable. Simple. Cost Effective
making your process easy, quick and affordable.

The T-Classic 9 is a 9-liter lightweight, portable steam autoclave designed as a cost-effective sterilization solution. Fully automatic with an easy one touch operation, delivering short cycles, starting at **18 Minutes!***

Ideally designed to meet portability needs of a small practitioner or clinic.

The T-Classic 9 can also be used in combination with another autoclave or as a back up.

T-Classic 9 allows you to have your instruments sterilized and ready quickly!

*Cycle times may vary, depending on load and weight.



Tuttnauer performance with the simplicity and affordability your small practice needs



Fast

The T-Classic 9 offers competitively quick unwrapped cycle times, in as short as 18 minutes*. The fast point-of-use sterilization makes your turnaround a priority, increasing efficiency, while eliminating bottlenecks.

*Cycle times may vary, depending on load & weight.



Intuitive & Simple

The T-Classic 9 is a simple one step operated device. With one touch, the T-Classic 9 will begin an automatic cycle until completion. An electronic detector will turn off the power if water level is insufficient. The perfect autoclave for limited space.



Compact & Portable

The T-Classic 9 is the lightest tabletop autoclave available, weighing only 11.7 lbs. allowing for easy portability, while requiring minimal counter space.



Safe

The interlock system prevents the lid from being opened while the vessel is pressurized. A thermal fuse is in place to protect against overheating and electrical faults.



Top view, easy access



Basket for easy loading and unloading



Easy to use control panel



Interlock safety feature

Technical Specifications

Model Number	T-Classic 9	Max load weight	3kg
Max Width / Height (overall)	13.4" / 13"	Sterilizing temperature	250°F
Max Net Weight	5.3 kg/11.7 lbs.	Min sterilizing time	18 min.
Chamber Diameter/Height	8.2"/9"	Cycle time - normal	28 min.
Chamber capacity	9 Liter / 2.4 Gallon	Voltage	120V
Max instrument length	9"	Watts	1200W
Max load weight	3kg/6.61 lbs.	Frequency	50-60Hz

Designed and manufactured to comply with the following standards and guidelines: ASME VIII Division 1, Medical Device Directive MDD 93/42/EEC, FDA 510(K) - Number K962903, CSA, UL, TSSA Canadian Design Registration

More products:



Ultrasonic Cleaners



Water Distiller



T-EDGE Autoclaves

LTR330-0121/T-ClassicUSA/REV-1/10.22