



PREVENTATIVE MAINTENANCE

**Protect Your Investment
Prevent Problems Before
They Start**

Tuttnauer Preventive Maintenance Kits are designed to keep your autoclave running smoothly year after year.

* Additional preventive maintenance kits are available beyond the options listed below.

MODEL#	PART #	COMPONENTS
NEW! T-TOP 10/11	PMTTOP 10/11	<ul style="list-style-type: none"> ▪ Hepa Filter ▪ System Clean (1 pk) ▪ Door Gasket ▪ PM Sticker
NEW! 3870HSG	PM3870HSG	<ul style="list-style-type: none"> ▪ Solenoid Valve (2ct.) ▪ Door Gasket ▪ Hepa Filter ▪ Steam Trap ▪ PM Sticker
NEW! Plazmax P50	PMPLZ50	<ul style="list-style-type: none"> ▪ Vacuum Pump Oil ▪ Vacuum Pump Filter ▪ Air Filter ▪ Oil Separator Filter ▪ P50 Gasket ▪ Peristaltic Hose for Pump ▪ H2O2 Vaporizer Filter ▪ Molecular Sieve for Odor Filter ▪ Plunger (2 ct.) ▪ O-Ring ▪ PM Sticker
NEW! Plazmax P110/160	PMPLZ110/160	<ul style="list-style-type: none"> ▪ Vacuum Pump Oil ▪ Vacuum Pump Filter ▪ Air Filter ▪ Oil Separator Filter ▪ P110/160 Gasket ▪ Peristaltic Hose for Pump ▪ H2O2 Vaporizer Filter ▪ Molecular Sieve for Odor Filter ▪ Plunger (2 ct.) ▪ O-Ring ▪ PM Sticker
NEW! 5596-1R	PM5596-1R	<ul style="list-style-type: none"> ▪ Steam Trap 1/2" W/O Bypass N125W ▪ Steam Trap 1/2" W/Bypass N125W ▪ Air Filter 0.01M 1/2" ▪ Door Gasket 9MM ▪ PM Sticker
NEW! 5596-1V	PM5596-1V	<ul style="list-style-type: none"> ▪ Steam Trap 1/2" W/O Bypass N125W ▪ Steam Trap 1/2" W/Bypass N125W ▪ Air Filter 0.01M 1/2" ▪ Door Gasket 20 X 20 ▪ PM Sticker
NEW! 6690-1V	PM6690-1V	<ul style="list-style-type: none"> ▪ Steam Trap 1/2" W/O Bypass N125W ▪ Steam Trap 1/2" W/Bypass N125W ▪ Air Filter 0.01M 1/2" ▪ Door Gasket 26 X 26 ▪ PM Sticker
NEW! 6690-2V	PM6690-2V	<ul style="list-style-type: none"> ▪ Steam Trap 1/2" W/O Bypass N125W ▪ Steam Trap 1/2" W/Bypass N125W ▪ Air Filter 0.01M 1/2" ▪ Door Gasket 26 X 26 (2ct.) ▪ PM Sticker
NEW! TIVA 8	119904671	<ul style="list-style-type: none"> ▪ Pre Filter AWD2 ▪ Impeller Washer (4 Ct) ▪ Flanged Bush (4 Ct) ▪ Hepa Filter ▪ Peristaltic Cartridge (2 Ct) ▪ PM Sticker
NEW! TIVA 2 Under Counter Or Above Counter	119904678	<ul style="list-style-type: none"> ▪ Pre Filter AWD2 ▪ Flanged Bush (4 Ct) ▪ Hepa Filter ▪ Peristaltic Cartridge (2 Ct) ▪ PM Sticker



PREVENTATIVE MAINTENANCE

**Protect Your Investment
Prevent Problems Before
They Start**

Tuttnauer Preventive Maintenance Kits are designed to keep your autoclave running smoothly year after year.

* Additional preventive maintenance kits are available beyond the options listed below.

MODEL #	PART #	COMPONENTS
2340M	PM2340M	<ul style="list-style-type: none"> ▪ Safety Valve ▪ Air Jet ▪ Door Gasket ▪ Chamber Brite (1 pk) ▪ PM Sticker
2540M/MK	PM2540M / PM2540MK	<ul style="list-style-type: none"> ▪ Safety Valve ▪ Air Jet ▪ Door Gasket ▪ Chamber Brite (1 pk) ▪ PM Sticker
3870EAP	PM3870EA	<ul style="list-style-type: none"> ▪ Safety Valve ▪ Air Jet ▪ Door Gasket ▪ Chamber Brite (1 pk) ▪ 6MM Plunger ▪ Hepa Filter ▪ 3MM Plunger (3 ct.) ▪ PM Sticker
3870M	PM3870M	<ul style="list-style-type: none"> ▪ Safety Valve ▪ Air Jet ▪ Door Gasket ▪ Chamber Brite (2 pks.) ▪ PM Sticker
1730M ValueKlave	PMVALUEKLAVE	<ul style="list-style-type: none"> ▪ Safety Valve ▪ Air Jet ▪ Door Gasket ▪ Chamber Brite (1 pk) ▪ PM Sticker