



# PREVENTATIVE MAINTENANCE

**Protect Your Investment  
Prevent Problems Before  
They Start**

Tuttnauer Preventive Maintenance Kits are designed to keep your autoclave running smoothly year after year.

\* Additional preventive maintenance kits are available beyond the options listed below.

MODEL#	PART #	COMPONENTS
<b>NEW!</b> T-TOP 10/11	PMTTOP 10/11	<ul style="list-style-type: none"> <li>Hepa Filter</li> <li>System Clean (1 pk)</li> <li>Door Gasket</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> 3870HSG	PM3870HSG	<ul style="list-style-type: none"> <li>Solenoid Valve (2ct.)</li> <li>Door Gasket</li> <li>Hepa Filter</li> <li>Steam Trap</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> Plazmax P50	PMPLZ50	<ul style="list-style-type: none"> <li>Vacuum Pump Oil</li> <li>Vacuum Pump Filter</li> <li>Air Filter</li> <li>Oil Separator Filter</li> <li>P50 Gasket</li> <li>Peristaltic Hose for Pump</li> <li>H2O2 Vaporizer Filter</li> <li>Molecular Sieve for Odor Filter</li> <li>Plunger (2 ct.)</li> <li>O-Ring</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> Plazmax P110/160	PMPLZ110/160	<ul style="list-style-type: none"> <li>Vacuum Pump Oil</li> <li>Vacuum Pump Filter</li> <li>Air Filter</li> <li>Oil Separator Filter</li> <li>P110/160 Gasket</li> <li>Peristaltic Hose for Pump</li> <li>H2O2 Vaporizer Filter</li> <li>Molecular Sieve for Odor Filter</li> <li>Plunger (2 ct.)</li> <li>O-Ring</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> 5596-1R	PM5596-1R	<ul style="list-style-type: none"> <li>Steam Trap 1/2" W/O Bypass N125W</li> <li>Steam Trap 1/2" W/Bypass N125W</li> <li>Air Filter 0.01M 1/2"</li> <li>Door Gasket 9MM</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> 5596-1V	PM5596-1V	<ul style="list-style-type: none"> <li>Steam Trap 1/2" W/O Bypass N125W</li> <li>Steam Trap 1/2" W/Bypass N125W</li> <li>Air Filter 0.01M 1/2"</li> <li>Door Gasket 20 X 20</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> 6690-1V	PM6690-1V	<ul style="list-style-type: none"> <li>Steam Trap 1/2" W/O Bypass N125W</li> <li>Steam Trap 1/2" W/Bypass N125W</li> <li>Air Filter 0.01M 1/2"</li> <li>Door Gasket 26 X 26</li> <li>PM Sticker</li> </ul>
<b>NEW!</b> 6690-2V	PM6690-2V	<ul style="list-style-type: none"> <li>Steam Trap 1/2" W/O Bypass N125W</li> <li>Steam Trap 1/2" W/Bypass N125W</li> <li>Air Filter 0.01M 1/2"</li> <li>Door Gasket 26 X 26 (2ct.)</li> <li>PM Sticker</li> </ul>
<b>NEW! TIVA 8</b>	119904671	<ul style="list-style-type: none"> <li>Pre Filter AWD2</li> <li>Impeller Washer (4 Ct)</li> <li>Flanged Bush (4 Ct)</li> <li>Hepa Filter</li> <li>Peristaltic Cartridge (2 Ct)</li> <li>PM Sticker</li> </ul>
<b>NEW! TIVA 2</b> Under Counter Or Above Counter	119904678	<ul style="list-style-type: none"> <li>Pre Filter AWD2</li> <li>Flanged Bush (4 Ct)</li> <li>Hepa Filter</li> <li>Peristaltic Cartridge (2 Ct)</li> <li>PM Sticker</li> </ul>
<b>NEW! TIVA 8</b> (3 Years Service)	119904674	<ul style="list-style-type: none"> <li>Pre Filter AWD2</li> <li>Impeller Washer (4 Ct)</li> <li>Chamber Gasket</li> <li>Suction Probe (orange)</li> <li>Suction Probe (blue)</li> <li>Flanged Bush (4 Ct)</li> <li>Hepa Filter</li> <li>Peristaltic Cartridge (2 Ct)</li> <li>PM Sticker</li> </ul>
<b>NEW! TIVA 2</b> Under Counter Only (3 Years Service)	119904683	<ul style="list-style-type: none"> <li>Chamber Gasket AWD</li> <li>Pre Filter Awd2</li> <li>Suction Probe (orange)</li> <li>Suction Probe (blue)</li> <li>Flanged Bush (4 Ct)</li> <li>Hepa Filter</li> <li>Peristaltic Cartridge (2 Ct)</li> <li>PM Sticker</li> </ul>
<b>NEW! TIVA 2</b> Counter Top Only (3 Years Service)	119904681	<ul style="list-style-type: none"> <li>Chamber Gasket AWD</li> <li>Pre Filter Awd2</li> <li>Flanged Bush (4 Ct)</li> <li>Hepa Filter</li> <li>Peristaltic Cartridge (2 Ct)</li> <li>PM Sticker</li> </ul>



# PREVENTATIVE MAINTENANCE

**Protect Your Investment  
Prevent Problems Before  
They Start**

Tuttnauer Preventive Maintenance Kits are designed to keep your autoclave running smoothly year after year.  
\* Additional preventive maintenance kits are available beyond the options listed below.

MODEL #	PART #	COMPONENTS
2340M	PM2340M	<ul style="list-style-type: none"> <li>▪ Safety Valve</li> <li>▪ Air Jet</li> <li>▪ Door Gasket</li> <li>▪ Chamber Brite (1 pk)</li> <li>▪ PM Sticker</li> </ul>
2540M/MK	PM2540M / PM2540MK	<ul style="list-style-type: none"> <li>▪ Safety Valve</li> <li>▪ Air Jet</li> <li>▪ Door Gasket</li> <li>▪ Chamber Brite (1 pk)</li> <li>▪ PM Sticker</li> </ul>
3870EAP	PM3870EA	<ul style="list-style-type: none"> <li>▪ Safety Valve</li> <li>▪ Air Jet</li> <li>▪ Door Gasket</li> <li>▪ Chamber Brite (1 pk)</li> <li>▪ 6MM Plunger</li> <li>▪ Hepa Filter</li> <li>▪ 3MM Plunger (3 ct.)</li> <li>▪ PM Sticker</li> </ul>
3870M	PM3870M	<ul style="list-style-type: none"> <li>▪ Safety Valve</li> <li>▪ Air Jet</li> <li>▪ Door Gasket</li> <li>▪ Chamber Brite (2 pks.)</li> <li>▪ PM Sticker</li> </ul>
1730M ValueKlave	PMVALUEKLAVE	<ul style="list-style-type: none"> <li>▪ Safety Valve</li> <li>▪ Air Jet</li> <li>▪ Door Gasket</li> <li>▪ Chamber Brite (1 pk)</li> <li>▪ PM Sticker</li> </ul>