

# Unveiling an innovative Virus protective shield

for our class B dental and  
medical autoclaves

Tuttnauer is proud to introduce a new pre-sterilization cycle to treat the viruses that may reside inside the autoclave chamber and provide best in class safety for the autoclave operators throughout the entire sterilization process.

For our larger autoclaves, we already implemented a filter to protect users from the air released into the environment. Now we are unveiling an innovative protective shield in our smaller autoclave models.



- Immediate delivery -



We have developed an innovative sterilization cycle that creates a protective shield which decontaminates the air within the autoclave, prior to its removal into the environment. This safe process is implemented through the entire sterilization cycle.



**This cycle has been developed based on academic research and literature, to ensure the complete and comprehensive elimination of viruses.**

*\*Available for the Elara 11 & 9 Tabletop Class B autoclaves.*

## Our fully automated sterilizer (Elara 11 & Elara 9) **Working for you**

### Clean Water & Waste Water Reservoirs

Reduces autoclave cleaning and maintenance and provides clean steam for every cycle

### Multi-Color Display

For easy reading, even from a distance

### Designed for Continuous Operation

Steam is continuously provided by an independent steam generator

### Complete Tractability


Archive cycle data including hard-copy printout

### Dynamic air removal technology

The high volume vacuum pump provides optimal sterilization results and removes moisture from the chamber. The pump also isolates the chamber from the outer environment during sterilization

### Closed Door Active Drying System

Maintains sterility and ensures efficient drying of packs and pouches



International Sales  
and Marketing  
E-mail: [info@tuttnauer.com](mailto:info@tuttnauer.com)  
[www.tuttnauer.com](http://www.tuttnauer.com)

Tuttnauer Europe b.v.  
Hoeksteen 11, 4815 PR  
PO Box 7191, 4800 GD Breda  
The Netherlands  
Tel: +31 765 423 510  
Fax: +31 765 423 540  
E-mail: [info@tuttnauer.nl](mailto:info@tuttnauer.nl)

Tuttnauer USA Co.  
25 Power Drive,  
Hauppauge, NY 11788  
Tel: +800 624 5836, +631 737 4850  
Fax: +631 737 0720  
E-mail: [info@tuttnauerUSA.com](mailto:info@tuttnauerUSA.com)

Distributed by:



Learn from our Experts  
Join our blog [tuttnauer.com/blog](http://tuttnauer.com/blog)

